

DID YOU KNOW ?

If you have Zoll or Physio-Control defibrillators that you can plug your Philips ECG patient cable into them by using our input adapter?



ME590030 allows the use of Philips Healthcare 12 pin patient cable with AAMI Standard 6 pin patient ECG input - pigtail

ECG Input Adapter. Allows the use of a Philips Healthcare monitor 12 pin patient ECG cable with an AAMI Standard 6 pin (3 lead or 5 lead) patient ECG ...

Our Price: \$187.00

[Add to Wish List](#) | [Add to Basket](#) | [More info...](#)



ME590030K-5 allows the use of Philips Healthcare 5 wire 12 pin patient ECG cables with the LIFEPAK 12,15 or 20

ECG Input Adapter for the Physio-Control LP12, LP15 or LP20. Allows the use of a Philips Healthcare 5 wire 12 pin patient ECG cable with the Physi...

Our Price: \$174.00

[Add to Wish List](#) | [Add to Basket](#) | [More info...](#)



ME590030K3 allows the use of Philips Healthcare 3 wire 12 pin patient ECG cables with the LIFEPAK 12, 15 or 20

ECG Input Adapter for the Physio-Control LP12, LP15 and LP20. Allows the use of a Philips Healthcare 3 wire 12 pin patient ECG cable with the Physi...

Our Price: \$194.00

[Add to Wish List](#) | [Add to Basket](#) | [More info...](#)



ME590030X allows the use of Philips Healthcare 12 Lead patient ECG cables with Propaq MD & Zoll X Series ECG input- pigtail

ECG Input Adapter. Allows the use of a Philips Healthcare 12 Pin - 12 Lead cable with a Propaq MD and Zoll X Series 12 lead patient ECG input. Thi...

Our Price: \$197.00

[Add to Wish List](#) | [Add to Basket](#) | [More info...](#)

MAGUIREENTERPRISES.COM

800-548-9686