

The 1/4" Phone Jack Outputs referred to in our cable descriptions connect to the devices listed below:

Manufacturer	Type	Model	ECG/BP	Special Notes
Arrow-Kontron	Balloon Pumps	Auto Cat2 Wave	ECG	
Arrow-Kontron	Balloon Pumps	Auto Cat2 Wave	BP	
Datascope/ Maquet	Balloon Pumps	CS300, Cardio Save	ECG	Combined outputs standard with ECG on tip, BP on ring, individual 1/4" outputs for ECG & IBP optional
Datascope/ Maquet	Balloon Pumps	CS300, Cardio Save	BP	Combined outputs standard with ECG on tip, BP on ring, individual 1/4" outputs for ECG & IBP optional
Datex	Monitors	F-CU8 processors with ECG Synchro Box	ECG	
Philips/HP	Monitors	All	ECG	
Philips/HP	Monitors	All	BP	1/8" output optional
Philips/HP	Defibrillators	HeartStart XL, HeartStart MRx, Codemaster XL, Codemaster XL+, Codemaster 100	ECG	
Physio-Control	Defibrillators	LP6, LP7, LP9	ECG	
Quinton	Hemodynamic	Q-Cath, SCS IV	ECG	3 circuit plug, BP on tip, ECG on ring
Quinton	Hemodynamic	Q-Cath, SCS IV	BP	3 circuit plug, BP on tip, ECG on ring
Siemens	Hemodynamic	Cathcor, Cathcor LX, Recor	ECG	
Siemens	Hemodynamic	Cathcor, Cathcor LX, Recor	BP	

MAGUIRE Enterprises, Inc.
1-800-548-9686

a member of The HackettGroup of Companies